

ScheBo® • Pancreas Elastase 1 Quick™ Canine



Use

The ScheBo® • Pancreas Elastase 1 Quick™ Canine is a visual immunochromatographic rapid test for detecting canine pancreatic elastase 1 in stool samples for non-invasive exocrine pancreatic function testing.

Main indications:

- Diagnosis/exclusion of pancreatic involvement in association with Gastrointestinal Symptoms, Maldigestion, Steatorrhea, Poor Hair Coat, Diarrhea, Weight Loss
- Diagnosis/exclusion of exocrine pancreatic insufficiency caused by e.g. Chronic Pancreatitis, Diabetes Mellitus, Pancreatic Cancer, Pancreatic Acinar Atrophy, etc.

Advantages:

- Based on monoclonal antibodies, substitution therapy has no influence on the test result
- Absolutely pancreas-specific
- No blood sampling
- No starvation periods needed
- Intestinal inflammation and blockage of the pancreatic ducts do not interfere with the results
- High specificity and sensitivity

Test Principle

Canine pancreatic elastase 1 (cPE 1) is a proteolytic digestive enzyme produced and secreted by the pancreas, which remains undegraded during intestinal transit and can be determined in stool. Its concentration in feces directly reflects the exocrine pancreatic function.

The ScheBo® • Pancreas Elastase 1 Quick™ Canine test is based on an immunochromatographic method. Canine pancreatic elastase 1 is detected by two specific monoclonal antibodies. Canine pancreatic elastase 1 in the stool sample reacts with a monoclonal antibody bound to gold particles. This complex migrates along the test membrane and reaches the test line (T), which has a second monoclonal antibody against pancreatic elastase 1 attached.

If the exocrine pancreatic function is normal (= high canine pancreatic elastase 1 concentration), the gold-labeled antibody + canine pancreatic elastase 1 complex binds to the test line (T) and a pink color becomes visible. In the event of an exocrine pancreatic insufficiency (= low canine elastase 1 concentration) the sample contains no antibody + canine elastase 1 complex that can bind to the test line (T) so no color becomes visible.

Development of a pink control line (C) guarantees that sample application and migration have taken place correctly and that the test was properly performed.

 **ScheBo • Biotech AG**

Netanyastrasse 3
35394 Giessen
Germany
Phone: +49-641-49960
www.schebo.com

Storage and Shelf Life

Test

The test must be stored at +4°C to +27°C and brought to room temperature just prior to use if necessary.

Stool sample

After taking the stool sample, it can be stored at room temperature for up to 5 days. Within these 5 days, either the test must be performed or the sample frozen at -20°C for longer-term storage. The deep-frozen sample is stable for up to 1 year.

Interferences

Very watery stools can lead to false results due to a dilution effect.

 **In-vitro diagnostic test**

 **Not for re-use**


- Do not use after the expiry date.
- Do not use reagents from different lots.
- Only open the test cassette's packaging shortly before performing the test.
- Make sure the foil bag containing the test cassette is undamaged.
- All samples should be regarded as potentially infectious. Therefore, wear disposable gloves when conducting the test and dispose of the samples, extracts and test cassettes accordingly.
- The extraction buffer contains small quantities of sodium azide.

Contents of the kit



- ① 1 instructions for use
- ② 1 ScheBo® • Pancreas Elastase 1 Quick™ Canine ready-to-use stool extraction system consisting of:
 - a) yellow dosing tip
 - b) yellow cone
 - c) pre-filled tube (extraction buffer with detergent and sodium azide (< 0.05%))
- ③ 1 Pipette
- ④ 1 ScheBo® • Pancreas Elastase 1 Quick™ Canine cassette (packed in an aluminum sachet)


Test Procedure



1. Turn the yellow dosing tip of the extraction system anti-clockwise and remove the yellow dosing tip by pulling it up.



2. Insert the yellow dosing tip in three different places in the stool sample to a depth of 1 cm (all notches must be filled with stool).




3. Please check that all notches of the yellow dosing tip are filled with stool.

Test Procedure



4. Insert the yellow dosing tip with the stool through the extraction system and turn the tip clockwise to close it.



5. Shake well and tap the tube, if necessary, until all the stool has been removed from the notches in the dosing tip.



6. Leave to stand for 10 minutes.



7. Give the tube a final shake. Caution: no stool should remain attached to the yellow dosing tip. If stool still remains stuck to the dosing tip, the extraction system can be left to stand for up to 1 hour in order to free the stool by repeated shaking.



8. Tear open the aluminum packaging and remove the test cassette.



9. Remove the stool sample extract from the extraction system with a pipette.



10. Using the pipette apply 4 drops of stool extract into the circular sample well on the cassette.



11. Wait exactly 10 minutes and then read the results. Results which are read later may be false.

Interpretation of the Test Results

Quality Control: The test contains an inbuilt control. When a pink band develops in the control region (C) the test has been performed correctly.

Normal Two pink bands develop, one in the control region (C) and one in the test region (T).

A high concentration of canine pancreatic elastase1 indicates a normal exocrine pancreatic function. The test stripe (T) must be clearly recognizable as a line, although it may be weaker than that of the control (C).



Low One pink band appears in the control region (C). No band develops in the test region (T). **A low concentration of canine pancreatic elastase1 indicates exocrine pancreatic insufficiency.**



Invalid If no pink band appears, the test has not worked.



Application and Test Interpretation

The test is used to measure canine pancreatic elastase1 in stool samples. A low concentration of canine pancreatic elastase1 indicates exocrine pancreatic insufficiency.

Reference concentration:

- normal canine pancreatic elastase1 concentration (= Values > 10 µg elastase1/g stool) is an indicator for a normal exocrine pancreatic function.
- low canine pancreatic elastase1 concentration (= Values < 10 µg elastase1/g stool) can be an indicator for an exocrine pancreatic insufficiency.

Two pink lines:

A normal level of canine pancreatic elastase1 was found in the stool sample at the time the test was performed. High values of canine pancreatic elastase1 concentrations are standard values and are an indicator for a normal exocrine pancreatic function. In case of disease symptoms, further diagnostic investigations should be conducted.

One pink line:

A low level of canine pancreatic elastase1 was found in the stool sample at the time the test was performed. Low canine pancreatic elastase1 concentrations indicate exocrine pancreatic insufficiency.

Treatment Option:

An exocrine pancreatic insufficiency can be treated with enzyme substitution therapy.

Performance:

The Pancreas Elastase1 Quick Canine test has > 95% sensitivity and specificity when compared to the Pancreatic Elastase1 Canine ELISA.

References

BATTERSBY I. A., PETERS, I. R., DAY M. J., WORMAN A. J., HALL E. J. (2005) Effect of intestinal inflammation on fecal elastase concentration in dogs *Vet Clin Pathol* 34 (1): 49-51

MORGAN J. A., MOORE L. E. (2009) A quick review of canine pancreatic insufficiency, *DVM* 360

SPILLMANN T. (1996) Zur Diagnostik der exokrinen Pankreasinsuffizienz beim Hund – Möglichkeiten und Grenzen der Labordiagnostik chronischer Pankreastopathien. *Kleintierpraxis* 41: 653-68

SPILLMANN T., CHAUDHRY, Y.S., WIBERG M., EIGENBRODT E. SZIGOLEIT A., WESTERMARCK E. (1999) Aktueller Stand der Labordiagnostik chronischer Pankreaserkrankungen beim Hund – Ein Methodenvergleich. Abstract von der 45. Jahrestagung der Deutschen Veterinärmedizinischen Gesellschaft Fachgruppe Kleintierkrankheiten, 07.- 10.10. 1999, Giessen

SPILLMANN T., EIGENBRODT E., SZIGOLEIT A. (1998) Die Bestimmung und klinische Relevanz der fäkalen pankreatischen Elastase beim Hund. *Tierärztl. Praxis* 26 (K): 364-368

SPILLMANN T., EIGENBRODT E., SZIGOLEIT A. (1999) Faecal Pancreatic Elastase in Dogs – Determination and Diagnostic Value. Abstract of the 7th Annual Congress of the European Society of Veterinary Internal Medicine, September 11 – 13, 1997, Lyon, France

SPILLMANN T., EIGENBRODT E., SZIGOLEIT A., CHAUDHRY Y.S. (1998) The clinical experiences with a species-specific ELISA for canine faecal pancreatic elastase. Abstract of the 44. Jahrestagung d. Deutschen Veterinärmedizinischen Gesellschaft Fachgruppe Kleintierkrankheiten, 19.-22.11.98 in Stuttgart, Germany and abstract of the Proceedings of the 2nd Annual Conference of the European Society of Veterinary Comparative Nutrition, Vienna, Austria

SPILLMANN T., MÜLLER E., GRÜNBAUM E.-G. (1994) Pancreatic Stimulation Test – A new technique for the diagnosis of exocrine pancreatic insufficiency. Focus on Gastroenterology Proceedings of the WALTHAM/ESCG Symposia: 55-7

SPILLMANN T., WIBERG M. E., TEIGELKAMP S., FAILING K., CHAUDHRY Y.S., KIRSCH A., EIFLER R., WESTERMARCK E., EIGENBRODT E., SZIGOLEIT A. (2000) Canine Faecal Pancreatic Elastase (cE1) in Dogs with Clinical Exocrine Pancreatic Insufficiency, Normal Dogs and Dogs with Chronic Enteropathies, *The European Journal of Comparative Gastroenterology* 5 (2): 1-6

SPILLMANN T., WIBERG M., TEIGELKAMP S., EIGENBRODT E., SZIGOLEIT A., CHAUDHRY Y.S., EIFLER R., WESTERMARCK E. (2000) Faecal Pancreatic Elastase (cE1) in Normal Dogs, Dogs with Chronic Enteritis and with Exocrine Pancreatic Insufficiency, Abstract of The Annual Congress of the European Society of Veterinary Internal Medicine (ESVIM), September 14-19, 2000, Neuchatel, Switzerland

SPILLMANN T., WITTKER A., TEIGELKAMP S., EIM C., BURKHARDT E., EIGENBRODT E., SZIGOLEIT A., (2001) An Immunoassay for Canine Pancreatic Elastase 1 as an Indicator for Exocrine Pancreatic Insufficiency in Dogs *J Vet Diagn Invest* 13:468-474

SPILLMANN T. (2003) Labordiagnostik von Pankreaserkrankungen des Hundes – Möglichkeiten und Grenzen. *Tierärztl Prax* 31 (K): 49-59

STROMBECK D.R. AND GUILDFORD G.W. (1993) *Small Animal Gastroenterology*, 2nd Edn. Wolfe Publishing Ltd., London 296-319, 344-356, 429-458