



SAFETY DATA SHEET

1. Product and Company Identification

Material name	Plasma, free hemoglobin Calibrator, Controls and Linearity Check Set.
Version #	2
Issue date	08/01/14
Revision date -	03/02/16
Supersedes date -	
CAS #	Mixture
Kit number	C464-10,11,12,0L
Product use	For the IN VITRO diagnostic use in calibration and quality control of serum assays
Manufacturer information	
Corporate Headquarters	Catachem 353 Christian Street Oxford, CT 06478 www.catacheminc.com Phone: +1203-262-0330
Emergency Telephone Numbers	+1203-232-1979

2. Hazards Identification

Physical state	Liquid.
Appearance	Clear to yellow liquid.
Emergency overview	Health injuries are not known or expected under normal use. Contains treated human serum base materials which should be treated as potentially hazardous and should be handled using good laboratory practices. The serum used has been tested for the Anti-HCV 3, ANTI-hiv-1/2, HBsAG and STS (RPR) by test kits approved by the FDA. The treated serum was found negative for HIV-1 and HCV by PCR assay.
OSHA regulatory status	This product is not hazardous according to OSHA 29CFR 1910.1200.
Potential health effects	
Routes of exposure	Skin contact. Eye contact.
Eyes	May cause eye irritation.
Skin	May cause skin irritation.
Inhalation	In high concentrations, vapors may be irritating to the respiratory system.
Ingestion	May cause discomfort if swallowed.
Target organs	Eyes.
Chronic effects	No data available.
Signs and symptoms	Direct contact with skin and eyes may cause irritation.
Potential environmental effects	The product components are not classified as environmentally hazardous. However, this does not exclude the possibility that large or frequent spills can have a harmful or damaging effect on the environment.

3. Composition / Information on Ingredients

The manufacturer lists no ingredients as hazardous according to OSHA 29 CFR 1910.1200.

4. First Aid Measures

First aid procedures

Eye contact In case of contact, immediately flush eyes with fresh water for at least 15 minutes while holding the eyelids open. Remove contact lenses if worn. Get medical attention if irritation persists.

Skin contact For skin contact flush with large amounts of water while removing contaminated clothing. Get medical attention if irritation develops and persists.

Inhalation Move to fresh air. Call a physician if symptoms develop or persist.

Ingestion If material is ingested, immediately contact a poison control center.

Notes to physician Provide general supportive measures and treat symptomatically.

General advice If you feel unwell, seek medical advice.

5. Fire Fighting Measures

Flammable properties This product is not flammable.

Extinguishing media

Suitable extinguishing media Extinguish with water spray, carbon dioxide, dry chemical or material appropriate for the surrounding fire.

Unsuitable extinguishing media None known.

Protection of firefighters

Specific hazards arising from the chemical None known.

Protective equipment and precautions for firefighters Self-contained breathing apparatus and full protective clothing must be worn in case of fire.

Fire fighting

equipment/instructions Use standard firefighting procedures and consider the hazards of other involved materials.

Hazardous combustion

Products None known.

6. Accidental Release Measures

Personal precautions Keep unnecessary personnel away. Wear appropriate protective equipment and clothing during clean-up. Do not touch damaged containers or spilled material unless wearing appropriate protective clothing.

Environmental precautions Do not allow to enter drains, sewers or watercourses.

Methods for containment Absorb spillage with non-combustible, absorbent material.

Methods for cleaning up Absorb spill with vermiculite or other inert material. Dispose of waste in accordance with all applicable federal, state, local and provincial environmental regulations, per Section 13.

Other information Clean up in accordance with all applicable regulations.

7. Handling and Storage

Handling	Avoid contact with skin and eyes. Wash thoroughly after handling. In case of insufficient ventilation, wear suitable respiratory equipment. Handle and open container with care.
Storage	Store at 2-8°C (35-46°F). Store in a closed container away from incompatible materials.

8. Exposure Controls / Personal Protection

Occupational exposure limits	No exposure limits noted for ingredient(s).
Exposure guidelines	Follow standard monitoring procedures.
Engineering controls	Provide adequate ventilation. Observe Occupational Exposure Limits and minimize the risk of inhalation of vapors. Provide easy access to water supply and eye wash facilities.
Personal protective equipment	
Eye / face protection	Wear approved safety glasses or goggles.
Skin protection	Wear lab coat or other protective garments. Remove contaminated clothing promptly.
Respiratory protection	Under normal conditions, respirator is not normally required.
General hygiene considerations	Handle in accordance with good industrial hygiene and safety practice.

9. Physical & Chemical Properties

Appearance	Colorless to yellow liquid.
Physical state	Liquid.
Form	Liquid.
Color	Colorless to yellow
Odor	Odorless.
Odor threshold	Not available.
pH	7.0-8.0 @ (25°C)
Vapor pressure	Not available.
Vapor density	Not available.
Boiling point	Not available.
Melting point/Freezing point	Not applicable.
Solubility (water)	Soluble.
Specific gravity	1-1.1
Flash point	Not available.
Flammability limits in air, upper, % by volume	Not available.
Flammability limits in air, lower, % by volume	Not available.
Auto-ignition temperature	Not available.
Evaporation rate	Not available.

10. Chemical Stability & Reactivity Information

Chemical stability	Material is stable under normal conditions.
Conditions to avoid	No data available.
Incompatible materials	No data available.
Hazardous decomposition	
Products	None known.
Possibility of hazardous	
Reactions	Hazardous polymerization does not occur.

11. Toxicological Information

Sensitization	Not classified.
Acute effects	May cause discomfort if swallowed.
Local effects	Ingestion may cause irritation and nausea.
Chronic effects	No data available.
Carcinogenicity	Not classifiable as to carcinogenicity to humans.
Epidemiology	No epidemiological data is available for this product.
Mutagenicity	Not classified.
Reproductive effects	Not classified.
Symptoms and target organs	Direct contact with skin and eyes may cause irritation.
Further information	No other specific acute or chronic health impact noted.

12. Ecological Information

Ecotoxicity	The product is not classified as environmentally hazardous. However, this does not exclude the possibility that large or frequent spills can have a harmful or damaging effect on the environment.
Environmental effects	An environmental hazard cannot be excluded in the event of unprofessional handling or disposal.
Persistence and degradability	No data is available on the degradability of this product.
Bioaccumulation / Accumulation	Not available.
Mobility in environmental media	The product is completely soluble in water.

13. Disposal Considerations

Disposal instructions	Contaminated instruments and surfaces should be disinfected in accordance with your employer's chemical-specific and universal/standard precautions.
Waste from residues / unused products	Dispose in accordance with all applicable regulations.
Contaminated packaging	Empty containers should be taken to an approved waste handling site for recycling or disposal.

14. Transport Information

DOT	Not regulated as a hazardous material by DOT.
------------	---

IATA Not regulated as dangerous goods.
IMDG Not regulated as dangerous goods.
TDG Not regulated as dangerous goods.

15. Regulatory Information

US federal regulations This product is not hazardous according to OSHA 29CFR 1910.1200. This mixture is a component of an in vitro diagnostic device regulated by the U.S. Food and Drug Administration.

TSCA Section 12(b) Export Notification (40 CFR 707, Subpart. D)

Not regulated.

Clean Air Act (CAA) Section 112 Hazardous Air Pollutants (HAPs) List

Not regulated.

CERCLA (Superfund) reportable quantity (lb.) (40 CFR 302.4)

None

Superfund Amendments and Reauthorization Act of 1986 (SARA)

Hazard categories Immediate Hazard - No
 Delayed Hazard - No
 Fire Hazard - No
 Pressure Hazard - No
 Reactivity Hazard - No

Section 302 extremely hazardous substance (40 CFR 355, Appendix A) Section 311/312 (40 CFR370) No

Drug Enforcement

Administration (DEA) (21 CFR

1308.11-15) Not controlled

Canadian regulations This product has been classified in accordance with the hazard criteria of the CPR and the MSDS contains all the information required by the CPR.

WHMIS status Non-controlled

Inventory status

Country(s) or region Inventory name On inventory (yes/no)*

Australia	Australian Inventory of Chemical Substances (AICS)	Yes
Canada	Domestic Substances List (DSL)	Yes
Canada	Non-Domestic Substances List (NDSL)	No
China	Inventory of Existing Chemical Substances in China (IECSC)	Yes
Europe	European Inventory of Existing Commercial Chemical Substances (EINECS)	Yes
Europe	European List of Notified Chemical Substances (ELINCS)	No
Japan	Inventory of Existing and New Chemical Substances (ENCS)	Yes
Korea	Existing Chemicals List (ECL)	Yes
New Zealand	New Zealand Inventory	Yes
Philippines	Philippine Inventory of Chemicals and Chemical Substances (PICCS)	

United States
& Puerto Rico Toxic Substances Control Act (TSCA) Inventory

Yes

*A "Yes" indicates this product complies with the inventory requirements administered by the governing country(s)

State regulations This product does not contain a chemical known to the State of California to cause cancer, birth defects or other reproductive harm.

US - California Proposition 65 - Carcinogens & Reproductive Toxicity (CRT): Listed substance

Not listed.

US. Massachusetts RTK - Substance List

Not regulated.

US. New Jersey Worker and Community Right-to-Know Act

Not regulated.

US. Pennsylvania RTK - Hazardous Substances

Not regulated.

Mexico regulations This safety data sheet was prepared in accordance with the Official Mexican Standard (NOM-018-STPS-2000).

16. Other Information

Recommended restrictions Use in accordance with supplier's recommendations.

Further information HMIS® is a registered trade and service mark of the NPCA.

HMIS® ratings

Health: 0

Flammability: 0

Physical hazard: 0

NFPA ratings

Health: 0

Flammability: 0

Instability: 0

Disclaimer The above information is believed to be correct, but does not purport to be all inclusive and should be used only as a guide. Catachem Inc. shall not be held liable for any damage resulting from handling or from contact with the above product.