

Small footprint.



Big potential.

The CatachemWell-T analyzer does what the “big cats” can with minimal hassle and a fraction of the cost.

Reliable, accurate chemistry tests at your fingertips.





CATACHEM

CHEMWELL-T



Improve Quality of Care:

- Easy to operate, calibrate twice a month
- Results within 15 minutes of startup
- Runs routine and most specialty tests
- Perfect STAT or back-up analyzer
- Complete QC package
- No proprietary consumables
- Compact footprint with robust design
- Patient database software included
- Easy LIS integration
- Customizable reporting styles
- User-friendly, programmable software
- Simple economical maintenance

Technical Features:

Throughput:

Up to 100 end point reactions per hour.

Up to 75 kinetic reactions per hour.

Reagents and Samples:

48 spaces for reagents and samples.

2 uL minimum sample volume.

Wavelengths:

Customizable based on lab's needs.

Thermal Control:

Cuvette carrier 37°C or room temperature.

Refrigerated reagent rack.

Cuvettes:

Reusable plastic.

8 strips of 5 cuvettes for 40 reactions.



CATACHEM